

VET-JECT®

18-LL Needles

- ▶ Copper content of hub: 59%
- ▶ Generously-shaped metal hubs, 18mm length, Luer Lock
- ▶ Hubs are marked with mm and gauge information
- ▶ ISO color-coded label for quick identification of diameter
- ▶ Packaging: 12pcs per case, 144pcs per box



Gauge	Dimensions	Color Code
23 G x 3/8"	0.60 x 10 mm	Blue
23 G x 5/8"	0.60 x 15 mm	Blue
23 G x 3/4"	0.60 x 20 mm	Blue
22 G x 3/8"	0.70 x 10 mm	Black
22 G x 5/8"	0.70 x 15 mm	Black
22 G x 3/4"	0.70 x 20 mm	Black
21 G x 3/8"	0.80 x 10 mm	Green
21 G x 5/8"	0.80 x 15 mm	Green
21 G x 1"	0.80 x 25 mm	Green
20 G x 3/8"	0.90 x 10 mm	Yellow
20 G x 1/2"	0.90 x 13 mm	Yellow
20 G x 5/8"	0.90 x 15 mm	Yellow
19 G x 3/8"	1.00 x 10 mm	Light Yellow
19 G x 1/2"	1.00 x 13 mm	Light Yellow
19 G x 5/8"	1.00 x 15 mm	Light Yellow
19 G x 3/4"	1.00 x 20 mm	Light Yellow
19 G x 1"	1.00 x 25 mm	Light Yellow
19 G x 1 1/4"	1.00 x 30 mm	Light Yellow
18 G x 3/8"	1.20 x 10 mm	Light Red
18 G x 1/2"	1.20 x 13 mm	Light Red
18 G x 5/8"	1.20 x 15 mm	Light Red
18 G x 3/4"	1.20 x 20 mm	Light Red
18 G x 1"	1.20 x 25 mm	Light Red
18 G x 1 1/4"	1.20 x 30 mm	Light Red
18 G x 1 1/2"	1.20 x 40 mm	Light Red
-	1.30 x 15 mm	Pink
-	1.30 x 20 mm	Pink
17 G x 3/8"	1.40 x 10 mm	Dark Red
17 G x 1/2"	1.40 x 13 mm	Dark Red
17 G x 5/8"	1.40 x 15 mm	Dark Red
17 G x 3/4"	1.40 x 20 mm	Dark Red
17 G x 1"	1.40 x 25 mm	Dark Red
17 G x 1 3/8"	1.40 x 35 mm	Dark Red
17 G x 1 1/2"	1.40 x 40 mm	Dark Red
-	1.50 x 15 mm	Dark Red
-	1.50 x 20 mm	Dark Red

Gauge	Dimensions	Color Code
-	1.50 x 25 mm	Dark Red
-	1.50 x 30 mm	Dark Red
-	1.50 x 35 mm	Dark Red
-	1.50 x 40 mm	Dark Red
-	1.50 x 50 mm	Dark Red
16 G x 3/8"	1.60 x 10 mm	White
16 G x 1/2"	1.60 x 13 mm	White
16 G x 5/8"	1.60 x 15 mm	White
16 G x 3/4"	1.60 x 20 mm	White
16 G x 1"	1.60 x 25 mm	White
16 G x 1 1/4"	1.60 x 30 mm	White
16 G x 1 3/8"	1.60 x 35 mm	White
16 G x 1 1/2"	1.60 x 40 mm	White
15 G x 1/2"	1.80 x 13 mm	Dark Blue
15 G x 5/8"	1.80 x 15 mm	Dark Blue
15 G x 3/4"	1.80 x 20 mm	Dark Blue
15 G x 1"	1.80 x 25 mm	Dark Blue
15 G x 1 1/4"	1.80 x 30 mm	Dark Blue
15 G x 1 3/8"	1.80 x 35 mm	Dark Blue
15 G x 1 1/2"	1.80 x 40 mm	Dark Blue
15 G x 1 3/4"	1.80 x 45 mm	Dark Blue
14 G x 1/2"	2.00 x 13 mm	Light Green
14 G x 5/8"	2.00 x 15 mm	Light Green
14 G x 3/4"	2.00 x 20 mm	Light Green
14 G x 1"	2.00 x 25 mm	Light Green
14 G x 1 1/4"	2.00 x 30 mm	Light Green
14 G x 1 3/8"	2.00 x 35 mm	Light Green
14 G x 1 1/2"	2.00 x 40 mm	Light Green
14 G x 1 3/4"	2.00 x 45 mm	Light Green
14 G x 2"	2.00 x 50 mm	Light Green
-	2.20 x 15 mm	Dark Green
-	2.20 x 20 mm	Dark Green
-	2.20 x 25 mm	Dark Green
-	2.20 x 30 mm	Dark Green
-	2.20 x 35 mm	Dark Green
13 G x 1 3/8"	2.40 x 35 mm	Pink

Other sizes are available according to customer's requests

