


BERPU-JECT®

Insulin Pen Needles

- ▶ Precision lancet point grinding, both tips are silicone treated for smooth penetration
- ▶ Thin wall technology ensures optimal insulin flow
- ▶ Exquisite design of screw helps the rotation more comfortable
- ▶ Various needle sizes are suitable for all types of patients
- ▶ Compatible for a wide range of insulin pens in the market



Code	Dimensions	Color Code
IN3404	0.18x4mm 34Gx4mm	
IN3304	0.20x4mm 33Gx4mm	
IN3204	0.23x4mm 32Gx4mm	
IN3205	0.23x5mm 32Gx5mm	
IN3206	0.23x6mm 32Gx6mm	
IN3104	0.25x4mm 31Gx4mm	
IN3105	0.25x5mm 31Gx5mm	
IN3106	0.25x6mm 31Gx6mm	
IN3108	0.25x8mm 31Gx8mm	
IN3005	0.30x5mm 30Gx5mm	
IN3008	0.30x8mm 30Gx8mm	
IN3010	0.30x10mm 30Gx10mm	
IN2906	0.33x6mm 29Gx6mm	
IN2910	0.33x10mm 29Gx10mm	
IN2912	0.33x12mm 29Gx12mm	



Insulin Pen	Brand
Novo Nordisk	InDuo®/InnoLet®/NovoPen®/NovoLet®/Innovo®/Flex Pen®
Eli Lilly	Humulin Pen®/Humalog Pen®
Sanofi-Aventis	Lantus Solostar Pen®/OptiPen®/OptiClik®
Owen Mumford	AutoPen®
Becton Dickinson	BD Pen®
B.Braun	Omnican Pen®
Berlin Chemie	BerliPen®
Haselmeier	DiaPen®