


BERPU-JECT®

Cannulas

- ▶ According to ISO 9626:2016 standard
- ▶ Made of medical grade stainless steel
- ▶ Fully automatic production line ensures high quality and precision
- ▶ Low cobalt or other materials are also available according to customer's requirements



Gauge	OD MIN	OD MAX	ID MIN (RW)	ID MIN (TW)
0.18 (34G)	0.178	0.191	0.064	0.091
0.20 (33G)	0.203	0.216	0.089	0.105
0.23 (32G)	0.229	0.241	0.089	0.105
0.25 (31G)	0.254	0.267	0.114	0.125
0.30 (30G)	0.298	0.320	0.133	0.165
0.33 (29G)	0.324	0.351	0.133	0.190
0.36 (28G)	0.349	0.370	0.133	0.190
0.40 (27G)	0.400	0.420	0.184	0.241
0.45 (26G)	0.440	0.470	0.232	0.292
0.50 (25G)	0.500	0.530	0.232	0.292
0.55 (24G)	0.550	0.580	0.280	0.343
0.60 (23G)	0.600	0.673	0.317	0.370
0.70 (22G)	0.698	0.730	0.390	0.440
0.80 (21G)	0.800	0.830	0.490	0.547
0.90 (20G)	0.860	0.920	0.560	0.635
1.10 (19G)	1.030	1.100	0.648	0.750
1.20 (18G)	1.200	1.300	0.790	0.910
1.40 (17G)	1.400	1.510	0.950	1.156
1.60 (16G)	1.600	1.690	1.100	1.283
1.80 (15G)	1.750	1.900	1.300	1.460
2.10 (14G)	1.950	2.150	1.500	1.600
2.40 (13G)	2.300	2.500	1.700	1.956
2.70 (12G)	2.650	2.850	1.950	2.235
3.00 (11G)	2.950	3.150	2.200	2.464
3.40 (10G)	3.300	3.500	2.500	2.819

Other sizes are available according to customer's requests

