

# VET-JECT®

## Veterinary IV Catheters

- ▶ PP tubing material
- ▶ Small and fixed wings increase patient comfort
- ▶ Thin-walled catheter with a large inner diameter assures good flow
- ▶ Popular size 26G is specially for small species, including reptiles and birds



Gauge	Outer Diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code
26G	0.60	15	12	
24G	0.70	19	22	
22G	0.90	25	36	
20G	1.00	32	58	
18G	1.20	32	90	
17G	1.50	45	120	
16G	1.80	45	172	
14G	2.00	45	208	

