

# BERPU-JECT<sup>®</sup>

## Blood Collection Needles

- ▶ Ultra sharp, siliconized triple bevel cannula minimizes pain during injection
- ▶ Excellent self-sealing rubber sleeve prevents blood flowing out
- ▶ Color-coded cap and hub assure accurate and convenient gauge identification
- ▶ Provide fast and convenient blood multiple sampling by single puncture



Code	Dimensions	Color Code
BN2325	0.60x25mm 23Gx1"	
BN2225	0.70x25mm 22Gx1"	
BN2238	0.70x38mm 22Gx1 1/2"	
BN2125	0.80x25mm 21Gx1"	
BN2138	0.80x38mm 21Gx1 1/2"	
BN2025	0.90x25mm 20Gx1"	
BN2038	0.90x38mm 20Gx1 1/2"	
BN1838	1.20x38mm 18Gx1 1/2"	



## Visible Blood Collection Needles

- ▶ Smart patent design offers more convenient operation and more effective puncture
- ▶ Self-made thin wall needle makes the blood collection more quickly and reduces the incidence rate of hemolysis
- ▶ The vein blood flow is visible right after puncture can solve the difficulty of collecting blood

Code	Dimensions	Color Code
VBN2225	0.70x25mm 22Gx1"	
VBN2238	0.70x38mm 22Gx1 1/2"	
VBN2125	0.80x25mm 21Gx1"	
VBN2138	0.80x38mm 21Gx1 1/2"	
VBN2025	0.90x25mm 20Gx1"	
VBN2038	0.90x38mm 20Gx1 1/2"	

