

BERPU-CARE[®]


Safety Needles

- ▶ According to ISO 7864 standard
- ▶ Safety mechanism is an integrated part of the needle
- ▶ Comfortable and easy operation, single-handed activation
- ▶ Bevel oriented needle for easy injection
- ▶ ISO color-coded hub for easy gauge identification



Code	Dimensions	Color Code
SN3012	0.30x12mm 30Gx1/2"	
SN2716	0.40x16mm 27Gx5/8"	
SN2620	0.45x20mm 26Gx3/4"	
SN2520	0.50x20mm 25Gx3/4"	
SN2425	0.55x25mm 24Gx1"	
SN2325	0.60x25mm 23Gx1"	
SN2330	0.60x30mm 23Gx1 1/4"	
SN2230	0.70x30mm 22Gx1 1/4"	
SN2140	0.80x40mm 21Gx1 1/2"	
SN2040	0.90x40mm 20Gx1 1/2"	
SN1940	1.10x40mm 19Gx1 1/2"	
SN1840	1.20x40mm 18Gx1 1/2"	

Other sizes are available according to customer's requests



Syringes With Safety Needles

- ▶ Comfortable and easy operation, single-handed activation
- ▶ Can prevent puncture accident effectively
- ▶ Simple combination structure is more economical than safety syringes

Syringe	Needle Gauge	Needle Length
1ml, 3ml, 5ml, 10ml, 20ml, 50ml Luer Slip / Luer Lock	18G~30G	13~50mm

