

# VET-JECT®

## Oral Syringes


- ▶ Deliver food or drug to animals
- ▶ Material: PP, silicone rubber, silicone oil
- ▶ O-Ring plunger design for smooth and accurate delivery
- ▶ Visible dosing marks allow easy measurement
- ▶ Non-Luer tip, can not fit with hypodermic needles



Code	Size
-	0.5ml, 1ml, 3ml, 5ml, 10ml, 20ml

## TPX Syringes

- ▶ Version A for easy dosing and version B for more precision
- ▶ Highly visible barrel scale with extended dose increments
- ▶ Exchangeable metal Luer Lock nozzle
- ▶ Fully boilable design for disinfection
- ▶ Individually packaged including spare o-ring parts
- ▶ Distinctive retail outer packaging
- ▶ Fitting drench nozzle available as accessory
- ▶ Dosage volume indicated on piston rod for Version B
- ▶ Dosage nut on the piston rod to lock in the dose for Version B



Code	Size
-	10ml, 20ml, 30ml, 50ml

## Veterinary Syringes

- ▶ Hard pack design, more convenient and safe, easily disposed and transportable
- ▶ Heat sealed cap for assurance of sterility

Code	Size
-	3ml, 5ml, 10ml, 20ml, 30ml, 60ml

