

VET-JECT®

21-LL/GW Needles

- ▶ Copper content of hub: 59%
- ▶ Distinctive bison bevel for long service life
- ▶ Hubs are available in two types: Luer Lock and Gewind (thread), 21mm length
- ▶ Thick-walled needle tube with exceptional stability and elasticity
- ▶ Generously-shaped metal hubs are marked with mm and gauge information
- ▶ Direct removal from transparent container ensures reduced risk of injury and contamination
- ▶ ISO color-coded label for quick identification of diameter
- ▶ Packaging: 10pcs per case, 100pcs per box



Gauge	Dimensions	Color Code
21 G x 3/8"	0.80 x 10 mm	Green
20 G x 3/8"	0.90 x 10 mm	Yellow
20 G x 5/8"	0.90 x 15 mm	Yellow
19 G x 3/8"	1.00 x 10 mm	Light Green
19 G x 5/8"	1.00 x 15 mm	Light Green
19 G x 3/4"	1.00 x 20 mm	Light Green
19 G x 1"	1.00 x 25 mm	Light Green
18 G x 3/8"	1.20 x 10 mm	Light Red
18 G x 5/8"	1.20 x 15 mm	Light Red
18 G x 3/4"	1.20 x 20 mm	Light Red
18 G x 1"	1.20 x 25 mm	Light Red
18 G x 1 3/8"	1.20 x 35 mm	Light Red
17 G x 3/8"	1.40 x 10 mm	Dark Red
17 G x 5/8"	1.40 x 15 mm	Dark Red
17 G x 3/4"	1.40 x 20 mm	Dark Red
17 G x 1"	1.40 x 25 mm	Dark Red
17 G x 1 1/4"	1.40 x 30 mm	Dark Red
17 G x 1 3/8"	1.40 x 35 mm	Dark Red
16 G x 3/8"	1.60 x 10 mm	White
16 G x 5/8"	1.60 x 15 mm	White
16 G x 3/4"	1.60 x 20 mm	White
16 G x 1"	1.60 x 25 mm	White

Gauge	Dimensions	Color Code
16 G x 1 1/4"	1.60 x 30 mm	White
16 G x 1 3/8"	1.60 x 35 mm	White
16 G x 1 1/2"	1.60 x 40 mm	White
15 G x 5/8"	1.80 x 15 mm	Dark Blue
15 G x 3/4"	1.80 x 20 mm	Dark Blue
15 G x 1"	1.80 x 25 mm	Dark Blue
15 G x 1 1/4"	1.80 x 30 mm	Dark Blue
15 G x 1 3/8"	1.80 x 35 mm	Dark Blue
15 G x 1 1/2"	1.80 x 40 mm	Dark Blue
14 G x 5/8"	2.00 x 15 mm	Light Green
14 G x 3/4"	2.00 x 20 mm	Light Green
14 G x 1"	2.00 x 25 mm	Light Green
14 G x 1 1/4"	2.00 x 30 mm	Light Green
14 G x 1 3/8"	2.00 x 35 mm	Light Green
14 G x 1 1/2"	2.00 x 40 mm	Light Green
14 G x 2"	2.00 x 50 mm	Light Green
-	2.20 x 25 mm	Dark Green
-	2.20 x 30 mm	Dark Green
-	2.20 x 35 mm	Dark Green
13G x 1"	2.40 x 25 mm	Pink
13G x 1 1/4"	2.40 x 30 mm	Pink

Other sizes are available according to customer's requests

