


VET-JECT®

18-LL Aqua Needles

- ▶ Copper content of hub: 59%
- ▶ Needle length tolerance is +/-0.25mm controlled by projector
- ▶ Generously-shaped metal hubs, 18mm length, Luer Lock
- ▶ Hubs are marked with mm and gauge information
- ▶ ISO color-coded label for quick identification of diameter
- ▶ Packaging: 12pcs per case, 144pcs per box





Gauge	Dimensions	Color Code
25G	0.50 x 3mm, 4mm, 5mm	
23G	0.60 x 3mm, 4mm, 5mm	
22G	0.70 x 3mm, 4mm, 10mm, 15mm	



Other sizes are available according to customer's requests

14-LL Needles

- ▶ Copper content of hub: 59%
- ▶ Hub length 14mm, Luer Lock
- ▶ Hubs are marked with mm and gauge information
- ▶ ISO color-coded label for quick identification of diameter
- ▶ Packaging: 12pcs per case, 144pcs per box



Gauge	Dimensions	Color Code
23 G x 3/8"	0.60 x 10 mm	
23 G x 5/8"	0.60 x 15 mm	
23 G x 3/4"	0.60 x 20 mm	
22 G x 3/8"	0.70 x 10 mm	
22 G x 5/8"	0.70 x 15 mm	
22 G x 3/4"	0.70 x 20 mm	
21 G x 3/8"	0.80 x 10 mm	
21 G x 5/8"	0.80 x 15 mm	
21 G x 1"	0.80 x 25 mm	
20 G x 3/8"	0.90 x 10 mm	
20 G x 1/2"	0.90 x 13 mm	
20 G x 5/8"	0.90 x 15 mm	
19 G x 3/8"	1.00 x 10 mm	
19 G x 1"	1.00 x 25 mm	
19 G x 1 1/4"	1.00 x 30 mm	
18 G x 1/2"	1.20 x 13 mm	
18 G x 1"	1.20 x 25 mm	
18 G x 1 1/4"	1.20 x 30 mm	
18 G x 1 1/2"	1.20 x 40 mm	

Gauge	Dimensions	Color Code
17 G x 3/8"	1.40 x 10 mm	
17 G x 1/2"	1.40 x 13 mm	
17 G x 3/4"	1.40 x 20 mm	
17 G x 1"	1.40 x 25 mm	
17 G x 1 1/2"	1.40 x 40 mm	
-	1.50 x 15 mm	
16 G x 3/8"	1.60 x 10 mm	
16 G x 1"	1.60 x 25 mm	
16 G x 1 1/2"	1.60 x 40 mm	
15 G x 1/2"	1.80 x 13 mm	
15 G x 1"	1.80 x 25 mm	
15 G x 1 3/8"	1.80 x 35 mm	
15 G x 1 1/2"	1.80 x 40 mm	
14 G x 1/2"	2.00 x 13 mm	
14 G x 1"	2.00 x 25 mm	
14 G x 1 3/8"	2.00 x 35 mm	
-	2.20 x 15 mm	
-	2.20 x 35 mm	
13 G x 1 3/8"	2.40 x 35 mm	

Other sizes are available according to customer's requests