

# BERPU-JECT®

## Blunt Fill Needles

- ▶ Designed for medication draw-up and re-constitution
- ▶ 45° bevel for easy penetration into vials and avoids needle stick injury
- ▶ Standard hub fits on all sizes of standard syringes


### Blunt Fill Needles With Filter

- ▶ 5micron filter membrane is designed to prevent larger particles being drawn into the syringe when preparing medications from ampoules
- ▶ 0.2micron filter membrane is designed to prevent air during using



Code	Dimensions	Color Code
BFN1840	18GX1 1/2"	
BFN18405	18GX1 1/2" 5micron	
BFN184002	18GX1 1/2" 0.2micron	
BFN1850	18Gx2"	
BFN18505	18Gx2" 5micron	
BFN185002	18Gx2" 0.2micron	

Other sizes are available according to customer's requests



## Dispensing Needles

- ▶ Self-made cannula for dispensing use
- ▶ Special designed needle tips can avoid generating fragments when inserted in a bottle stopper



Code	Dimensions	Color Code
BDN1830	1.20x30mm 18Gx1 1/4"	
BDN1840	1.20x40mm 18Gx1 1/2"	
BDN1735	1.50x35mm 17Gx1 2/5"	
BDN1730	1.50x30mm 17Gx1 1/4"	
BDN1750	1.50x50mm 17Gx2"	
BDN1625	1.60x25mm 16Gx1"	
BDN1630	1.60x30mm 16Gx1 1/4"	
BDN1640	1.60x40mm 16Gx1 1/2"	

