

# VET-JECT®

## 14-LL Detectable Needles

- ▶ Copper content of hub: 59%
- ▶ Specifically designed with detectability and strength for meat industry
- ▶ A detectable alloy in the cannula material allows the needles to be as much as 100% detectable in metal detection systems
- ▶ Ultra sharp tri-bevelled tip for easier puncture, allows multiple injections without burring
- ▶ Hub length 14mm, Luer Lock
- ▶ Gauge stamped on hub for easy identification
- ▶ Single cartridge package, 100pcs per box



Size	Dimensions	Color Code
14D x 1"	2.10 x 25 mm	-
14D x 1 1/2"	2.10 x 40 mm	-
16D x 1/2"	1.82 x 13 mm	-
16D x 5/8"	1.82 x 15 mm	-
16D x 3/4"	1.82 x 20 mm	-
16D x 1"	1.82 x 25 mm	-
16D x 1 1/2"	1.82 x 40 mm	-
18D x 1/2"	1.47 x 13 mm	-
18D x 5/8"	1.47 x 15 mm	-
18D x 3/4"	1.47 x 20 mm	-
18D x 1"	1.47 x 25 mm	-
18D x 1 1/2"	1.47 x 40 mm	-
20D x 1/2"	1.06 x 13 mm	-
20D x 1"	1.06 x 25 mm	-
20D x 1 1/2"	1.06 x 40 mm	-

Other sizes are available according to customer's requests



## 14-LL Cartridge Needles

- ▶ Copper content of hub: 59%
- ▶ Hub length 14mm, Luer Lock
- ▶ Hubs are marked with mm and gauge information
- ▶ Single cartridge package, 100pcs per box



Gauge	Dimensions	Length
13G~23G	0.60~2.40mm	10~45mm

