

VET-JECT®

Aluminum Hub Needles

- ▶ Stainless steel cannula, ultra sharp, triple bevel
- ▶ Regular wall thickness is less likely to bend
- ▶ Aluminum hub, Luer Lock
- ▶ Single cartridge package, sterile
- ▶ Cartridges are ISO color-coded for easy gauge identification



Gauge	Dimensions	Length (mm)
20G	0.90mm	13 ≤ L ≤ 50
18G	1.20mm	
16G	1.60mm	
14G	2.00mm	



PP Hub With Aluminum Insert Needles

- ▶ Stainless steel cannula, ultra sharp, tri-bevelled
- ▶ Regular wall thickness is less likely to bend
- ▶ Metal insert increases needle strength for animal health applications
- ▶ PP hub is ISO color-coded for easy gauge identification
- ▶ Single blister package and cartridge package are available



Gauge	Dimensions	Length (mm)
30G	0.30mm	13 ≤ L ≤ 50
27G	0.40mm	
26G	0.45mm	
25G	0.50mm	
24G	0.55mm	
23G	0.60mm	
22G	0.70mm	
21G	0.80mm	
20G	0.90mm	
19G	1.10mm	
18G	1.20mm	

